



CYTOSTOR™

**PRESERVATION OF CELLS & VIRUSES
IN SUSPENSION**

CS-SNS; CS-LNS



Atelerix Handbook Series

Version CS01;1.1.1

TABLE OF CONTENTS

1. Components	1
1.1. Kit Contents	1
1.2. Components to be Supplied by the User.....	2
1.3. Before You Begin Using Cytostor™	2
2. Protocol overview	2
3. CytoStor™ (CS-SNS) for the Encapsulation of Cells	3
3.1. Gelation.....	3
3.2. Cell Density and Loading Guide	4
3.2.1. Cell Loading Guide.....	4
3.2.2. Recommended Cell Load for Different Cell Types	4
3.3. Release	5
3.4. Shipping your cells	5
4. CytoStor™ Large (CS-LNS) for the Encapsulation of Cells.....	6
4.1. Gelation.....	6
4.2. Cell Density and Loading Guide	7
4.2.1. Cell Loading Guide.....	7
4.2.2. Recommended Cell Load for Different Cell Types	7
4.3. Release	8
4.4. Shipping your cells	8
5. CytoStor™ (CS-SNS) for Encapsulation of Viruses.....	9
5.1. Gelation.....	9
5.2. Viral Titre and Loading Guide	10
5.2.1. Loading Guide	10
5.2.2. Recommended Viral Load	10
5.3. Release	11
5.4. Shipping your virus	11
6. Troubleshooting Guide	12
7. Statements	13
7.1. Kit Storage and Stability	13
7.2. Cellular Material	13
7.3. Trademarks	13
8. Appendix: Cell and Virus Storage Temperature Guide	13

1. COMPONENTS

1.1. KIT CONTENTS

PRODUCT CODE	COMPONENTS	UNITS	UNIT VOLUME	CELL SUSPENSION TO ADD
CS-SNS-03	CytoStor Vial (Beads)	3 Tubes	0.4ml	-
	Gel A (5x)	3 Tubes	0.12ml	0.48ml
	Dissolution Buffer	3 Tubes	1.1ml	-
CS-SNS-06	CytoStor Vial (Beads)	6 Tubes	0.4ml	-
	Gel A (5x)	6 Tubes	0.12ml	0.48ml
	Dissolution Buffer	6 Tubes	1.1ml	-
CS-SNS-12	CytoStor Vial (Beads)	12 Tubes	0.4ml	-
	Gel A (5x)	12 Tubes	0.12ml	0.48ml
	Dissolution Buffer	12 Tubes	1.1ml	-
CS-SNS-24	CytoStor Vial (Beads)	24 Tubes	0.4ml	-
	Gel A (5x)	24 Tubes	0.12ml	0.48ml
	Dissolution Buffer	24 Tubes	1.1ml	-
CS-SNS-50	CytoStor Vial (Beads)	50 Tubes	0.4ml	-
	Gel A (5x)	50 Tubes	0.12ml	0.48ml
	Dissolution Buffer	50 Tubes	1.1ml	-
CS-LNS-03	CytoStor Vial (Beads)	3 Tubes	2ml	-
	Gel A (5x)	3 Tubes	0.6ml	2.4ml
	Dissolution Buffer	3 Tubes	5ml	-
CS-LNS-06	CytoStor Vial (Beads)	6 Tubes	2ml	-
	Gel A (5x)	6 Tubes	0.6ml	2.4ml
	Dissolution Buffer	6 Tubes	5ml	-
CS-LNS-12	CytoStor Vial (Beads)	12 Tubes	2ml	-
	Gel A (5x)	12 Tubes	0.6ml	2.4ml
	Dissolution Buffer	12 Tubes	5ml	-
CS-LNS-24	CytoStor Vial (Beads)	24 Tubes	2ml	-
	Gel A (5x)	24 Tubes	0.6ml	2.4ml
	Dissolution Buffer	24 Tubes	5ml	-
CS-LNS-50	CytoStor Vial (Beads)	50 Tubes	2ml	-
	Gel A (5x)	50 Tubes	0.6ml	2.4ml
	Dissolution Buffer	50 Tubes	5ml	-

Note: Remove components from 2-8°C for at least 20 minutes before use.

1.2. COMPONENTS TO BE SUPPLIED BY THE USER

- 1000µl Pipettes and Tips
- Cell Culture Medium
- **Cells or viruses!!**

1.3. BEFORE YOU BEGIN USING CYTOSTOR™

1. Ensure CytoStor™ kits have not passed the expiry date stated on the packaging. Atelerix does not recommend using kits after this date.
2. Read the troubleshooting guide on page 12 to see our list of frequently asked questions. For any further queries, please email us at Sales@atelerix.co.uk.
3. Consult the Cell Density and Loading guide (page 4 or 7 for cells, page 10 for viruses) and the Cell and Virus Storage Temperature Guide on page 13.
4. CytoStor™ is intended for use solely in accordance with this protocol using the components provided within the kit.

2. PROTOCOL OVERVIEW



3. CYTOSTOR™ (CS-SNS) FOR THE ENCAPSULATION OF CELLS

3.1. GELATION

1. Ensure that all components are allowed to equilibrate to room temperature before use and that gels are at the bottom of their tubes. Conduct all steps in a laminar flow hood at room temperature.
2. Resuspend cells in the appropriate volume of culture medium (See 3.2 Cell Density and Loading Guide on page 4 for more information), ensuring thorough distribution of cells.
3. Add 0.48 mL of the cell culture suspension to the vial containing 0.12 mL of Gel A.
4. Gently mix until homogenous, with a pipette, ensuring that no bubbles are introduced (see troubleshooting guide on page 12).
5. Add 0.6 mL of the cells / Gel A mix to the CytoStor Vial¹.
6. Place the cap back on the tube and gently invert the gel / bead mixture several times until the beads are evenly distributed throughout the gel. Gently flick the tube to settle the contents, ensuring a tight seal (the gel will cure in situ within approximately 30 minutes, sample is ready to ship after 1 hour).
7. Store away from light in a polystyrene box at the recommended temperature for the cell type encapsulated. See the table on page 13 or, for the most up to date recommendations on storage temperatures and times, visit our [Compatibility section](#) on our website or contact sales@atelerix.co.uk.

For directions on shipping your samples, see **Section 3.4**.

¹Use the CytoStor Vial containing beads provided for encapsulation, storage, and release.



3.2. CELL DENSITY AND LOADING GUIDE

3.2.1. CELL LOADING GUIDE

Desired Cell Load Per Sample (x10 ⁶)	Cell Suspension Concentration (x10 ⁶ cells/mL)	Volume of Cell Suspension per Sample (mL)
1	2.08	0.48
2	4.17	0.48
3	6.25	0.48
4	8.33	0.48
5	10.42	0.48
6	12.50	0.48
7	14.58	0.48
8	16.67	0.48
9	18.75	0.48
10	20.83	0.48

For multiple encapsulations, multiply the volume of cell suspension (C) by the number of samples required.

3.2.2. RECOMMENDED CELL LOAD FOR DIFFERENT CELL TYPES

Cell Diameter	Example Cell Types	Recommended Cell Load per Encapsulation
4 – 10 Microns	Lymphocytes e.g T-Cells, PBMC's	2 x 10 ⁷ Cells
11 – 15 Microns	Fibroblasts CHO Cells HEK Cells	≤ 1 x 10 ⁷ Cells
16 – 30 Microns	Mesenchymal Stem Cells Monocytes Hepatocytes	≤ 1 x 10 ⁷ Cells

CytoStor™ has no lower limit to the number of cells which can be stored, however we recommend using enough cells to be easily recovered after centrifugation.

3.3. RELEASE

1. Ensure that all components and samples are allowed to equilibrate to room temperature before use and conduct all steps in a laminar flow hood at room temperature.
2. Using a pipette tip or syringe and needle, infuse 1 mL Dissolution Buffer into the bottom of the gel by piercing the gel.
3. Place the cap back on the tube and allow the gel to dissolve by occasionally agitating the tube by gentle inversion or rocking for 10 minutes (see troubleshooting guide on page 12).
4. Sediment cells by centrifugation at your cells' usual centrifugation settings (e.g. 350RCF for 5 minutes), remove supernatant, and re-suspend cells in medium of choice.

3.4. SHIPPING YOUR CELLS

Use appropriate controlled temperature packaging² when preparing cells for shipping to reduce the effect of ambient temperature change on the encapsulated cells during transit.

²For best cell recovery upon arrival, we recommend using the ICECATCH Solid Ambient or Cool Shipping Boxes.

Find out more & shop at <https://www.atelerix.co.uk/pages/variants-collection-page-accessories>



4. CYTOSTOR™ LARGE (CS-LNS) FOR THE ENCAPSULATION OF CELLS

4.1. GELATION

1. Ensure that all components are allowed to equilibrate to room temperature before use and that gels are at the bottom of their tubes. Conduct all steps in a laminar flow hood at room temperature.
2. Resuspend cells in the appropriate volume of culture medium (See 4.2 Cell Density and Loading Guide on page 7 for more information), ensuring thorough distribution of cells.
3. Add 2.4 mL of the cell culture suspension to the vial containing 0.6 mL of Gel A.
4. Gently mix until homogenous, with a pipette, ensuring that no bubbles are introduced (see troubleshooting guide on page 12).
5. Add 2.4 mL of the cells / Gel A mix to the CytoStor Vial¹.
6. Place the cap back on the tube and gently invert the gel / bead mixture several times until the beads are evenly distributed throughout the gel. Gently flick the tube to settle the contents, ensuring a tight seal (the gel will cure in situ within approximately 30 minutes, sample is ready to ship after 1 hour).
8. Store away from light in a polystyrene box at the recommended temperature for the cell type encapsulated. See the table on page 13 or, for the most up to date recommendations on storage temperatures and times, visit our [Compatibility section](#) on our website or contact sales@atelerix.co.uk.

For directions on shipping your samples, see **Section 4.4**.

¹Use the CytoStor Vial containing beads provided for encapsulation, storage, and release.

4.2. CELL DENSITY AND LOADING GUIDE

4.2.1. CELL LOADING GUIDE

Desired Cell Load Per Sample ($\times 10^6$)	Cell Suspension Concentration ($\times 10^6$ cells/mL)	Volume of Cell Suspension per Sample (mL)
5	2.08	2.40
10	4.17	2.40
15	6.25	2.40
20	8.33	2.40
25	10.42	2.40
30	12.50	2.40
35	14.58	2.40
40	16.67	2.40
45	18.75	2.40
50	20.83	2.40

For multiple encapsulations, multiply the volume of cell suspension (C) by the number of samples required.

4.2.2. RECOMMENDED CELL LOAD FOR DIFFERENT CELL TYPES

Cell Diameter	Example Cell Types	Recommended Cell Load per Encapsulation
4 – 10 Microns	Lymphocytes e.g T-Cells, PBMC's	1×10^8 Cells
11 – 15 Microns	Fibroblasts CHO Cells HEK Cells	$\leq 5 \times 10^7$ Cells
16 – 30 Microns	Mesenchymal Stem Cells Monocytes Hepatocytes	$\leq 5 \times 10^7$ Cells

CytoStor™ has no lower limit to the number of cells which can be stored, however we recommend using enough cells to be easily recovered after centrifugation.



4.3. RELEASE

1. Ensure that all components and samples are allowed to equilibrate to room temperature before use and conduct all steps in a laminar flow hood at room temperature.
2. Using a pipette or syringe and needle, infuse 5 mL Dissolution Buffer into the bottom of the gel by piercing the gel.
3. Place the cap back on the tube and allow the gel to dissolve by occasionally agitating the tube by gentle inversion or rocking for 10 minutes (see troubleshooting guide on page 12).
5. Sediment cells by centrifugation at your cells' usual centrifugation settings (e.g. 350RCF for 5 minutes), remove supernatant, and re-suspend cells in medium of choice.

4.4. SHIPPING YOUR CELLS

Use appropriate controlled temperature packaging² when preparing cells for shipping to reduce the effect of ambient temperature change on the encapsulated cells during transit.

²For best cell recovery upon arrival, we recommend using the ICECATCH Solid Ambient or Cool Shipping Boxes.

Find out more & shop at <https://www.atelrix.co.uk/pages/variants-collection-page-accessories>

5. CYTOSTOR™ (CS-SNS) FOR ENCAPSULATION OF VIRUSES

5.1. GELATION

1. Ensure that all components are allowed to equilibrate to room temperature before use and that gels are at the bottom of their tubes. Conduct all steps in a laminar flow hood at room temperature.
2. Resuspend virus in the appropriate volume of medium (See 5.2 Viral Titre and Loading Guide on page 10 for more information), ensuring thorough distribution of virus.
3. Add 0.48 mL of the virus suspension to the vial containing 0.12 mL of **Gel A**.
4. Gently mix until homogenous, with a pipette, ensuring that no bubbles are introduced (see troubleshooting guide on page 12).
5. Add 0.6 mL of the virus / **Gel A** mix to the **CytoStor Vial**¹.
6. Place the cap back on the tube and gently invert the **gel / bead** mixture several times until the beads are evenly distributed throughout the gel. Gently flick the tube to settle the contents, ensuring a tight seal (the gel will cure in situ within approximately 30 minutes, sample is ready to ship after 1 hour).
7. Store away from light in a polystyrene box at room temperature (15°C - 25°C).

For directions on shipping your samples, see **Section 5.4**.

¹Use the **CytoStor Vial** containing beads provided for encapsulation, storage, and release.



5.2. VIRAL TITRE AND LOADING GUIDE

5.2.1. LOADING GUIDE

Desired Viral Load per Sample (x 10 ⁸)	Virus Titre (x 10 ⁸ IFU/mL)	Volume of Virus Suspension per Sample (mL)
1	2.08	0.48
2	4.17	0.48
3	6.25	0.48
4	8.33	0.48
5	10.42	0.48
6	12.50	0.48
7	14.58	0.48
8	16.67	0.48
9	18.75	0.48
10	20.83	0.48

For multiple encapsulations, multiply the volume of viral suspension (C) by the number of samples required.

5.2.2. RECOMMENDED VIRAL LOAD

Sample	Example Types	Recommended Viral Load per Encapsulation
Virus	Lentivirus, Coronavirus, Viral Vectors	x 10 ⁹ IFU

CytoStor™ has no lower limit to the number of virus which can be stored, however we recommend using enough cells to be easily recovered after centrifugation.

5.3. RELEASE

1. Ensure that all components and samples are allowed to equilibrate to room temperature before use and conduct all steps in a laminar flow hood at room temperature.
2. Using a pipette tip or syringe and needle, infuse 1 mL **Dissolution Buffer** into the bottom of the gel by piercing the gel.
3. Place the cap back on the tube and allow the gel to dissolve by occasionally agitating the tube by gentle inversion or rocking for 10 minutes (see troubleshooting guide on page 12).
4. Virus can then either be:
 - a. Sedimented by ultracentrifugation. Following ultracentrifugation, remove supernatant and re-suspended virus in medium of choice.
 - b. Diluted in a medium of choice (minimum 1 in 4 dilution recommended).

5.4. SHIPPING YOUR VIRUS

Use appropriate controlled room temperature packaging² when preparing virus for shipping to reduce the effect of ambient temperature change on the encapsulated virus during transit.

²For best cell recovery upon arrival, we recommend using the ICECATCH Solid Ambient or Cool Shipping Boxes.

Find out more & shop at <https://www.atelierix.co.uk/pages/variants-collection-page-accessories>



6. TROUBLESHOOTING GUIDE

Problem / Question	Guidance
I have air bubbles in the gel after mixing with my media, is this a problem?	Air trapped within the gel will affect preservation, so bubbles should be eliminated before mixing with the beads. Allow time for the mixture to settle and the bubbles to travel to the surface before addition.
Can I use the kit to encapsulate organoids and / or spheroids?	We would recommend TissueReady™ for the encapsulation and storage of organoids and spheroids.
Can I ship the Dissolution Buffer in the same package as the samples?	Yes, the Dissolution Buffer is stable at a wide range of temperatures and can be shipped together with the encapsulated samples.
What are the recommended storage times and temperatures for my cell type?	A guide to the recommended storage times and temperatures can be found on our Compatibility section on our website. If you cannot find any recommendations for your cell type, please contact Sales@atelerix.co.uk .
Can I reuse the contents of the kit if I don't use it all?	No, there should only be sufficient volume for a set number of encapsulations per kit. Any spare reagents will not be sufficient to perform any additional encapsulations properly.
Can I use PBS instead of media when encapsulating samples?	No, PBS should not be used at any point as it inhibits and slowly reverses gelation.
Can I split the kit into smaller tubes to get more encapsulations?	No, we do not recommend removing the GelBase Beads from the tubes supplied or deviating from the volumes stated.
What if the beads have not fully dissolved after 10 minutes? Can the cells be allowed to sit in the buffer for longer?	This sometimes occurs when the medium used to store the beads fortifies the gel. Extend the incubation until the beads dissolve - this should not take much longer than 15 minutes. Please contact Sales@atelerix.co.uk for information on how long cells may be allowed to incubate in the Dissolution Buffer.

7. STATEMENTS

7.1. KIT STORAGE AND STABILITY

This kit is stable at 2-8°C for up to 12 months. Atelerix does not recommend using the kit after the expiry date stated on the packaging.

7.2. CELLULAR MATERIAL

Please ensure that cell cultures are free of fungal and bacterial contamination before proceeding.

7.3. TRADEMARKS

CytoStor™ is a trademark of Atelerix Ltd.

8. APPENDIX: CELL AND VIRUS STORAGE TEMPERATURE GUIDE

Cell Type	Recommended Storage Conditions	Tested Storage Time
PBMCs	2-8°C	3-Days
Jurkat Cells	15-25°C	14-Days
MSCs	15-20°C	21-Days
Monocytes	2-8°C	5-Days
Fibroblasts	15-20°C	14-Days
Viruses	15-25°C	14-Days

If you cannot find any recommendations for your cell type, please visit our [Compatibility](#) section on our website or contact Sales@atelerix.co.uk.



NOTES

